

InSpec™ Gamme Générale

DESIGNATION	CODE PRODUIT	CODE COULEUR	COMPOSITION	Concentration	NATURE	Bag-in-bottle	Filter in Trigger	Screw Cap	TISSUS	Sterile	Non Sterile	Gamma Irradiation	Aseptic Fill * See Note	STOCKAGE	In use shelf life
IPA 1L PREMIUM	IPAWFI1-1LS		Isopropyl Alcohol (Propan-2-ol)	70/30 v/v Alcohol to WFI Quality Water	Biocide	●				●		●		2Yrs	3mths
IPA 1L STANDARD	IPAWFI1-1LSTND				Biocide		●			●		●		2Yrs	3mths
IPA 500ML PREMIUM	IPAWFI-500MLS				Biocide	●				●		●		2Yrs	3mths
IPA 500ML STANDARD	IPAWFI1-500MLSTND				Biocide		●			●		●		2Yrs	3mths
IPA 5L	IPAWFI1-5LS				Biocide			●		●		●		2Yrs	10 hrs
IPA 1L NON-STERILE	IPAWFI-1LNS				Biocide		●			●		NA	NA	N/A	N/A
IPA 5L NON-STERILE	IPAWFI1-5LNS				Biocide			●		●		NA	NA	N/A	N/A
IPA 500ML NON-STERILE	IPAWFI-500MLNS				Biocide		●			●		NA	NA	N/A	N/A
IPA PREMIUM TUB WIPES	IPAWFI12TB130S				Biocide				●	●		●		2Yrs	N/A
IPA STANDARD TUB WIPES	IPAWFI12TB200S				Biocide				●	●		●		2Yrs	N/A
IPA STANDARD 15 POUCH WIPES	IPAWFIPCH15S				Biocide				●	●		●		2Yrs	10 hrs
IPA STANDARD 50 POUCH WIPES	IPAWFIPCH50S				Biocide				●	●		●		2Yrs	10 hrs
DE 1L PREMIUM	DEWFI10-1LS		Denatured Ethanol (also known as Industrial Methylated Spirits or IMS)	70/30 v/v Alcohol to WFI Quality Water	Biocide	●				●		●		2Yrs	3mths
DE 1L STANDARD	DEWFISTND-1LSTND				Biocide		●			●		●		2Yrs	3mths
DE 500ML PREMIUM	DEWFI11-500MLS				Biocide	●				●		●		2Yrs	3mths
DE 500ML STANDARD	DEWFI1-500MLSTND				Biocide		●			●		●		2Yrs	3mths
DE 5L	DEWFI10-5LS				Biocide			●		●		●		2Yrs	10 hrs
DE PREMIUM TUB WIPES	DEWFI12TB130S				Biocide				●	●		●		2Yrs	N/A
DE STANDARD TUB WIPES	DEWFI12TB200S				Biocide				●	●		●		2Yrs	N/A
DE STANDARD 15 POUCH WIPES	DEWFIPCH15S				Biocide				●	●		●		2Yrs	10 hrs
DE STANDARD 50 POUCH WIPES	DEWFIPCH50S				Biocide				●	●		●		2Yrs	10 hrs
OX 1L PREMIUM	OXPW40-1LS					Sterile Sporicidal Preparation	6% Hydrogen Peroxide 0.4% Peracetic Acid v/v	Sporicide	●				●		
OX 5L	OXPW40-5LS	Sporicide						●		●			●	2Yrs	10 hrs
OX BURSTABLE POUCH	OXPW40BP15S	Sporicide							●	●			●	2Yrs	10 hrs
QT 900MLS PREMIUM (**See Note)	QTWFI30-900MLS		Sterile Quaternary Ammonium Compound Disinfectant Cleaner	4.0% QUAT disinfectant to WFI	Biocide	●				●		●		2Yrs	3mths
QT 5L	QTWFI30-5LS				Biocide			●		●		●		2Yrs	10 hrs
QT BURSTABLE POUCH	QTWFI30BP15S				Biocide				●	●		●		2Yrs	2 hrs
QT CONCENTRATE	QTCNC30-100			Concentrate	Biocide			●		●		●		2Yrs	10 hrs
AN 900MLS PREMIUM (** See Note)	ANWFI20-900MLS		Sterile Diamine Disinfectant Cleaner	4.0% Diamine disinfectant to WFI	Biocide	●				●		●		2Yrs	3mths
AN 5L	ANWFI20-5LS				Biocide			●		●		●		2Yrs	10 hrs
AN BURSTABLE POUCH	ANWFI20BP15S				Biocide				●	●		●		2Yrs	2 hrs
AN CONCENTRATE	ANCNC20-100			Concentrate	Biocide			●		●		●		2Yrs	10 hrs
HA 1L PREMIUM	HAPW50-1LS		Sterile Sporicidal Preparation	Hypochlorous Acid 200-800ppm	Sporicide	●				●			●	1Yr	3mths
HA 5L PREMIUM	HAPW50-5LS				Sporicide			●		●			●	1Yr	10 hrs
HA BURSTABLE POUCH	HAPW50BP15S				Sporicide				●	●			●	1Yr	2 hrs
N10 1L MANUAL DETERGENT NON-STERILE	N10WFI1-1LS		Manual Detergent	4.0% Detergent to WFI	Detergent		●				●	NA	NA	N/A	N/A
N10 1L MANUAL DETERGENT STERILE	N10WFI1CC-1LS				Detergent	●				●		●		2Yrs	3mths
N10 5L MANUAL DETERGENT STERILE	N10WFI1-5LS				Detergent			●		●		●		2Yrs	10 hrs
N10 100ML CONCENTRATE	N10WFI1-100MLS			Concentrate	Detergent			●		●		●		2Yrs	10 hrs
N10 1L CONCENTRATE	N101-1CONC				Detergent			●		●		●		2Yrs	10 hrs

* Filled using Aseptic techniques – All packaging gamma irradiated at 25 to 35 kGy’s prior to fill in ISO 5 conditions
** 1 Litre bottle containing 900mls of disinfectant solution